





HiSorb Starter Kit (DVB/CWR/PDMS, short, CVD inert)

Part number: C-HSPKIT-04-D

Part number	Description	Quantity	Image (not to scale)
H4-DXABC-06	HiSorb probe, DVB/CWR/PDMS, CVD inert, short version, pk 6	1 (6 probes)	
C-HSPCCS	Cap with HiSorb septum, for 10 mL/20 mL crimp- top vials, magnetic, pk 100	1 (100 caps/septa)	
C-VCC20	Headspace vial, 20 mL, crimp-top, clear, round- bottomed, pk 100	1 (100 vials)	
C-HSPH	HiSorb handle, for manual probe extraction	1	
C-HSPVLP-6	HiSorb septum plug, pk 6	1 (6 plugs)	