

Centri Booster Pack

Part number: U-CENKIT

Part number	Description	Quantity	Image (not to scale)
C-DL010	DiffLok cap, stainless steel, 1/4", pk 10	6 (60 caps)	
U-COV07	O-ring, Size 007, pk 10	1 (10 O-rings)	
U-COV10	O-ring, Size 010, pk 10	1 (10 O-rings)	

Note

When using inert coated stainless-steel TD tubes, Markes recommend inert coated stainless steel DiffLok caps for the sampling end of the tube during automated analysis.

Associated part numbers:

- Inert coated stainless-steel DiffLok caps, pk 10 (Part number: C-DLS10)
- Inert coated stainless-steel DiffLok caps, pk 100 (Part number: C-DL1S0)
- Pair of DiffLok caps: one stainless-steel and one inert coated stainless-steel, pk 10 pairs (Part number C-DLP10)
- Pair of DiffLok caps: one stainless-steel and one inert coated stainless-steel, pk 100 pairs (Part number C-DL1P0)