

Automated TD Booster Pack (DiffLok cap pairs)

Part number: U-ULKIT-P

Part number	Description	Quantity	Image (not to scale)
C-DL010	DiffLok cap, stainless steel, 1/4", pk 10	8 (80 caps)	
C-DLS10	Cap, DiffLok, 1/4", Inert, pk 10	8 (80 caps)	
U-COV07	O-ring, Size 007, pk 10	1 (10 O-rings)	
U-COV10	O-ring, Size 010, pk 10	1 (10 O-rings)	

NOTE Markes recommends inert-coated stainless steel DiffLok caps for the sampling end of the tube during automated analysis.

Associated part numbers:

- Inert coated stainless-steel DiffLok caps, pk 10 (Part number: C-DLS10)
- Inert coated stainless-steel DiffLok caps, pk 100 (Part number: C-DL1S0)
- Pair of DiffLok caps: one stainless-steel and one inert coated stainless-steel, pk 10 pairs (Part number C-DLP10)
- Pair of DiffLok caps: one stainless-steel and one inert coated stainless-steel, pk 100 pairs (Part number C-DL1P0)